

# Choosing the proper gas and gas equipment in the laboratory, Part 2: Relationship between gas purity and laboratory analyzers

Michel Gastiger and Benjamin Jurcik

**G**as purity requirements are directly related to analyzer requirements and the field of use or the working domain. For example, Part 1 of this paper<sup>1</sup> illustrated that high-purity gases such as ALPHAGAZ 2 (AIR LIQUIDE Corp., Paris, France) and dedicated equipment must be used and implemented to analyze low-level concentration samples (e.g., <0.1 ppm) no matter which GC detectors are used. While it is evident that extremely high-purity gases can be used for measurements requiring a higher degree of precision, the reverse is not true. As a result, the authors' selection criteria for gas and material is to use the grades of gas and material necessary to maintain the good working condition of the analyzer and to provide the needed analytical precision.

This paper applies the same logic for the other most common analyzers in the laboratory. For each, the concentration domain to be analyzed will be indicated and the grade of gas recommended will be given. In regard to the grade of gas, dedicated equipment will be chosen.

Analyzers requiring gases have been segmented into the following areas: 1) emission spectroscopy, 2) absorption spectroscopy, 3) mass spectroscopy, 4) thermal analysis, 5) liquid chromatography or supercritical fluid chromatography, 6) surface analysis, 7) gas analysis, and 8) other analysis. Since gas chromatography is the largest analytical technique that uses gas, the relationship between gas purity and chromatography detectors has been described in detail.<sup>1</sup>

For each of the analytical techniques mentioned above, a grade of gas is recommended in relationship to the level of concentration of the sample.

A word of caution is needed: The levels of concentration (in mol/mol or wt/wt) of the sample to be analyzed are indications, not absolute values; there are always specific conditions that may arise requiring the use of a higher gas purity and that cannot be foreseen.

## Pure gases

In order to simplify the choice for the analyst, only two grades for helium, nitrogen, argon, hydrogen, oxygen, acetylene, nitrous oxide, and carbon dioxide (ALPHAGAZ 1 and ALPHAGAZ 2, AIR LIQUIDE Corp.) are proposed. The specifications of both pure gases are given in Table 1. The purity of these gases may be as high as 99.9995% (N<sub>2</sub>) or 99.99999% (H<sub>2</sub>). The large range of flow rate and pressure permits the gas to be fed to one or more analyzers. The typical specifications of

Table 1

Gases	Specifications of ALPHAGAZ 1 and ALPHAGAZ 2	
	ALPHAGAZ 1	ALPHAGAZ 2
Ar, H <sub>2</sub> , He, N <sub>2</sub>	H <sub>2</sub> O < 3 ppm O <sub>2</sub> < 2 ppm C <sub>n</sub> H <sub>m</sub> < 0.5 ppm	H <sub>2</sub> O < 0.5 ppm O <sub>2</sub> , C <sub>n</sub> H <sub>m</sub> , CO, CO <sub>2</sub> < 0.1 ppm N <sub>2</sub> < 0.1 ppm in Ar, He H <sub>2</sub> < 0.1 ppm in Ar, He, N <sub>2</sub>
Air	H <sub>2</sub> O < 3 ppm C <sub>n</sub> H <sub>m</sub> < 0.5 ppm	— H <sub>2</sub> O < 0.5 ppm
O <sub>2</sub>	H <sub>2</sub> O < 3 ppm C <sub>n</sub> H <sub>m</sub> < 0.5 ppm	H <sub>2</sub> , C <sub>n</sub> H <sub>m</sub> , CO, CO <sub>2</sub> < 0.1 ppm N <sub>2</sub> < 4 ppm NO <sub>x</sub> < 15 ppb (typical analysis) H <sub>2</sub> O < 3 ppm
CO <sub>2</sub>	H <sub>2</sub> O < 20 ppm	C <sub>n</sub> H <sub>m</sub> < 2 ppm O <sub>2</sub> < 2 ppm N <sub>2</sub> < 8 ppm
C <sub>2</sub> H <sub>2</sub>	N <sub>2</sub> < 0.4%	—
N <sub>2</sub> O	N <sub>2</sub> < 0.4%	—

Table 2

Gas	Purity	Specifications of ALPHAGAZ FLO generators		
		Maximum output flow rate (L/min)	Maximum output pressure (bar)	Typical analysis
Nitrogen	>99.99%	1–11	Up to 8.6	CO, CO <sub>2</sub> , O <sub>2</sub> < 1 ppm H <sub>2</sub> O < 2 ppm Hydrocarbons < 0.1 ppm
Hydrogen	>99.99999%	0.15–0.55	Up to 4	Not applicable
Air	Not applicable	1–30	Up to 8	Hydrocarbon (as CH <sub>4</sub> ) < 0.1 ppm

Table 3

ALPHAGAZ MIX	Range of ALPHAGAZ MIX	
	Concentrations	Uses
Methane–argon	5 or 10% CH <sub>4</sub> in argon	Electron capture detector, X-ray fluorescence, nuclear counter
Hydrogen–helium	40% H <sub>2</sub> in helium	Flame ionization detector
Hydrogen–argon	0.75–7% H <sub>2</sub> in argon	Spark emission

ALPHAGAZ FLO generators (AIR LIQUIDE Corp.) are given in Table 2. Analyzers that require special mixtures for working (e.g., X-ray fluorescence) are fed with ALPHAGAZ MIX grade (AIR LIQUIDE Corp.) (Table 3).

These specifications do not cover all the possible cases of critical impurities for a given analyzer or analysis; they are designed to ensure the quality of gas necessary for laboratory applications.

### Recommendations

The analyzers presented in this paper are primarily dedicated to the analysis of liquid or solid samples. The impurities that can be contained in gases are only present in gaseous phase and as a result do not interfere with the majority of analyses that will be done. For this reason, use of ALPHAGAZ 1 or ALPHAGAZ FLO is recommended for most analyzers.

Table 4

		Emission spectroscopy Detection level (mol/mol or wt/wt)			
Techniques	Gas %	<1000 ppm	<100 ppm	<10 ppm	<1 ppm
<b>Inductively coupled plasma (ICP)</b>					
Auxiliary plasma flow	Ar	ALPHAGAZ 1			
Plasma flow	Ar	ALPHAGAZ 1			
<b>Optical detector (ICP-OES)*</b>					
Purge of optical part	N <sub>2</sub>	ALPHAGAZ 1			
Purge of optical part	N <sub>2</sub>	ALPHAGAZ FLO			
<b>Mass spectrometer (ICP-MS)</b>					
Auxiliary plasma flow	Ar	N/A	ALPHAGAZ 1		
Plasma flow	Ar	N/A	ALPHAGAZ 1		
<b>Other</b>					
Organic solvent analysis	O <sub>2</sub>	ALPHAGAZ 1			
	H <sub>2</sub>	ALPHAGAZ 1			
<b>Spark emission</b>					
Plasma	Ar	ALPHAGAZ 1**	N/A		
Plasma	Ar + 1-7% H <sub>2</sub>	ALPHAGAZ MIX	N/A		
Purge of optical part	N <sub>2</sub>	ALPHAGAZ 1	N/A		
Purge of optical part	N <sub>2</sub>	ALPHAGAZ FLO	N/A		
<b>GD (glow discharge) OES</b>					
Plasma	Ar	N/A	ALPHAGAZ 2		
Purge of optical part	N <sub>2</sub>	N/A	ALPHAGAZ 1		
Purge of optical part	N <sub>2</sub>	N/A	ALPHAGAZ FLO		
<b>X-ray fluorescence</b>					
Flow counter	Ar-CH <sub>4</sub> (90-10)	ALPHAGAZ MIX			
Flow counter	C <sub>n</sub> H <sub>m</sub> mixtures in He	ALPHAGAZ 1			
Detector cooling	N <sub>2</sub>	Liquid nitrogen			
Liquid analysis	He	ALPHAGAZ 1			
<b>UV fluorescence</b>					
Purge of optical part	N <sub>2</sub>	N/A	ALPHAGAZ 1***		
Purge of optical part	N <sub>2</sub>	N/A	ALPHAGAZ FLO		
Analysis of Hg	Ar	N/A	ALPHAGAZ 1		
Analysis of H <sub>2</sub> S or sulfurs	Air	N/A	ALPHAGAZ FLO		
Analysis of H <sub>2</sub> S or sulfurs	Air	N/A	ALPHAGAZ 1*		
<b>Chemiluminescence or CDL (Chemiluminescence detector)</b>					
Process gas	O <sub>2</sub>	ALPHAGAZ 1			
Process gas	Air	ALPHAGAZ FLO			
Process gas	Air	ALPHAGAZ 1***			
O <sub>3</sub> analysis	Ethylene	N/A	N35		
<b>Flame spectrometry</b>					
Sample introduction	A <sub>2</sub>	ALPHAGAZ 1	N/A		
Flame	O <sub>2</sub>	ALPHAGAZ 1	N/A		
Flame	N <sub>2</sub> O	ALPHAGAZ 1	N/A		
Flame	C <sub>2</sub> H <sub>2</sub>	ALPHAGAZ 1	N/A		

\*OES, optical emission spectroscopy.

\*\*ALPHAGAZ 2 in the case of nitrogen analysis.

\*\*\*Special grade (e.g., POL or VEM) in the case of analysis below 1 ppm.

Table 5

Techniques	Gas	Absorption spectroscopy Detection level (mol/mol or wt/wt)			
		<1000 ppm	<100 ppm	<10 ppm	<1 ppm
<b>Atomic absorption spectroscopy (AAS) with graphite furnace (GFAAS)</b>					
Graphite furnace	Ar	ALPHAGAZ 1			
<b>Atomic absorption with flame (AAS or FAAS)</b>					
Flame C <sub>2</sub> H <sub>2</sub> /Air	Air	ALPHAGAZ 1			
Flame C <sub>2</sub> H <sub>2</sub> /Air	Air	ALPHAGAZ FLO			
Flame	C <sub>2</sub> H <sub>2</sub>	ALPHAGAZ 1			
Flame C <sub>2</sub> H <sub>2</sub> /N <sub>2</sub> O	N <sub>2</sub> O	ALPHAGAZ 1			
Flame C <sub>2</sub> H <sub>2</sub> /N <sub>2</sub> O	C <sub>2</sub> H <sub>2</sub>	ALPHAGAZ 1			
Hydride analysis	N <sub>2</sub>	ALPHAGAZ 1			
Hydride analysis	N <sub>2</sub>	ALPHAGAZ FLO			
<b>Infrared (IR) spectroscopy</b>					
Purge or zero gas	N <sub>2</sub>	ALPHAGAZ 1*			
Purge or zero gas	N <sub>2</sub>	ALPHAGAZ FLO			
<b>FTIR spectroscopy</b>					
Purge or zero gas	N <sub>2</sub>	ALPHAGAZ 1*			
Purge or zero gas	N <sub>2</sub>	ALPHAGAZ FLO			
Cooling of the MCT detector**	N <sub>2</sub>	Liquid nitrogen			
<b>Elemental analysis of C, O, N, H, and S</b>					
N analysis: TCD**	He, Ar	ALPHAGAZ 1	ALPHAGAZ 2		
O analysis: IR analyzer	He, Ar	ALPHAGAZ 1	ALPHAGAZ 2		
C analysis: IR analyzer	O <sub>2</sub>	ALPHAGAZ 1	ALPHAGAZ 2		
S analysis: IR analyzer	O <sub>2</sub>	ALPHAGAZ 1	ALPHAGAZ 2		
H analysis: TCD	Ar, N <sub>2</sub>	ALPHAGAZ 1	ALPHAGAZ 2		
<b>Photoacoustic spectroscopy (PAS)</b>					
N <sub>2</sub> for calibration		ALPHAGAZ 1			
	N/A				
<b>Nuclear magnetic resonance (NMR), electron paramagnetic resonance (EPR) or (ESR)</b>					
Cooling	He	Liquid helium			
Cooling	N <sub>2</sub>	Liquid nitrogen			
Cooling	Air	ALPHAGAZ FLO			
Sample spinning					
<b>RAMAN spectroscopy</b>					
Purge	N <sub>2</sub>	ALPHAGAZ 1			
Purge	N <sub>2</sub>	ALPHAGAZ FLO			
Purge	N <sub>2</sub>	Liquid nitrogen			
Cooling of the Ge detector					

\*ALPHAGAZ 2 in the case of CO<sub>2</sub> analysis (below 10 ppm).

\*\*MCT, mercury cadmium telluride; TCD, thermal conductivity detector.

Table 6

Techniques	Mass spectrometry	
	Gas	Grade
<b>MS</b>		
Fast atomic bombardment	Ar, Xe	ALPHAGAZ 1 or standard grades
Chemical ionization	NH <sub>3</sub> , CH <sub>4</sub> , isobutane	ALPHAGAZ 1 or standard grades
Tandem mass spectrometer (MS-MS)	Ar, N <sub>2</sub>	ALPHAGAZ 1 or standard grades
Atmospheric pressure ionization (API)	Ar, N <sub>2</sub>	ALPHAGAZ 1 or standard grades
<b>Coupling LC to MS</b>		
Venturi effect	Air, N <sub>2</sub>	ALPHAGAZ 1
Venturi effect	Air, N <sub>2</sub>	ALPHAGAZ FLO
Curtain gas	N <sub>2</sub> , He	ALPHAGAZ 1
Curtain gas	N <sub>2</sub>	ALPHAGAZ FLO
Nebulization	Air, N <sub>2</sub>	ALPHAGAZ 1
Nebulization	Air, N <sub>2</sub>	ALPHAGAZ FLO
<b>Coupling GC to MS</b>		
Open split coupling	He	ALPHAGAZ 1*

\*ALPHAGAZ 2 for low concentration analysis (below 10 ppm).

Table 7

Thermal analysis	Thermal analysis	
	Gas	Grade
<b>Thermogravimetric analysis (TGA), differential scanning, calorimetry (DSC)</b>		
Inert gases	Ar, N <sub>2</sub> , He	ALPHAGAZ 1
Process gases	O <sub>2</sub> , Air, H <sub>2</sub> , SO <sub>2</sub>	ALPHAGAZ 1 or standard grades
Process or inert gases	Air, N <sub>2</sub> , or H <sub>2</sub>	ALPHAGAZ FLO
Cooling	N <sub>2</sub>	Liquid nitrogen

Table 8

Techniques	Gas	Liquid or supercritical chromatography			
		Detection level (sample in solution)			
		g/L	g/L	g/L	mg/L
		0.1	0.01	10	<1
<b>Liquid chromatography (HPLC or LC)</b>					
Purge of solvent	He	ALPHAGAZ 1			
<b>Supercritical fluid chromatography (SFC)</b>					
Mobile phase	CO <sub>2</sub>	ALPHAGAZ 1 or SFE or SFC grades			
<b>Supercritical fluid extraction (SFE)</b>					
Mobile phase	CO <sub>2</sub>	ALPHAGAZ 1 or SFE or SFC grades			

Table 9

Surface analysis	Surface analysis	
	Gas	Grade
<b>X-ray photoelectron spectroscopy (XPS or electron spectroscopy for chemical analysis [ESCA])</b>		
Sputtering gas	Ar	ALPHAGAZ 1
Cooling of the detector	N <sub>2</sub>	Liquid nitrogen
<b>Auger electron spectroscopy (Auger or AES)*</b>		
Purging	N <sub>2</sub>	ALPHAGAZ 1
Cooling of the detector	N <sub>2</sub>	Liquid nitrogen
<b>Electron microscopy (SEM, TEM)*</b>		
Purging	N <sub>2</sub>	ALPHAGAZ 1
Cooling of EDAX detector	N <sub>2</sub>	Liquid nitrogen
<b>Specific surface (BET)*</b>		
Cooling	N <sub>2</sub>	Liquid nitrogen
Surface measurement	N <sub>2</sub>	ALPHAGAZ 2
Surface measurement	Kr	N48
Zero adjustment	He	ALPHAGAZ 2

\*AES, Auger electron spectroscopy; SEM, scanning electron microscopy; TEM, transmission electron microscopy; BET, Brunnauer, Emmett, and Teller method.

Nevertheless, for the analysis of very dilute samples, ALPHAGAZ 2 should be used to improve the baseline, particularly in cases involving the analysis of O, N, C in steel, or for the analysis of nitrogen in solids by glow discharge spectroscopy. For the analysis of atmospheric pollutants such as NO, NO<sub>2</sub>, or SO<sub>2</sub> by chemiluminescence or UV fluorescence, use of ALPHAGAZ 2 grade is suggested. For more assurance, certain grades such as pollution (POL) or Vehicle Emission Zero (VEM), which include specifications on these pollutants, may also be used. Similarly, two special grades of CO<sub>2</sub> have been developed for supercritical fluid chromatography applications (SFE [supercritical fluid extraction] and SFC [supercritical fluid chromatography] grades) that contain very low levels of heavy hydrocarbons. Except in these particular cases, ALPHAGAZ 1 and ALPHAGAZ FLO are usually the grades recommended.

Tables 4–11 list the most common analyzers using gas and the recommended grade. Only two grades of gas are necessary for supplying the most common laboratory analyzers. Furthermore, ALPHAGAZ FLO provides an alternative solution for supplying the gases, which are well adapted for some analyzers.

Table 10

Techniques	Gas	Other analysis			
		Detection level (sample in solution)			
		g/L	g/L	g/L	mg/L
		0.1	0.01	10	<1
<b>Water analysis (TOC, TOX)*</b>					
Oxidation or flushing	O <sub>2</sub>	ALPHAGAZ 1		ALPHAGAZ 2	
Purging and desorption	N <sub>2</sub> , He, Ar	ALPHAGAZ 1		ALPHAGAZ 2	
Purging and desorption	N <sub>2</sub>	ALPHAGAZ FLO			
<b>Measurement of water in solids (Karl Fischer analysis)</b>					
Water desorption	Air, N <sub>2</sub>	ALPHAGAZ 1		ALPHAGAZ 2	
Water desorption	Air, N <sub>2</sub>	ALPHAGAZ FLO			
<b>Metal analysis in petrochemistry (Wickbold method)</b>					
Flame	O <sub>2</sub>	ALPHAGAZ 1			
Flame	H <sub>2</sub>	ALPHAGAZ 1			
Flame	H <sub>2</sub>	ALPHAGAZ FLO			
Stripping	N <sub>2</sub>	ALPHAGAZ 1			
Stripping	N <sub>2</sub>	ALPHAGAZ FLO			

\*TOC, total organic carbon; TOX, total organic halides.

Table 11

Techniques	Gas	%	Gas analysis			
			Detection level (mol/mol)			
			<1000 ppm	<100 ppm	<10 ppm	<1 ppm
<b>Hygrometry</b>						
Purging or zero gas	N <sub>2</sub>	N/A	ALPHAGAZ 1		ALPHAGAZ 2	
Purging or zero gas	N <sub>2</sub>	N/A	ALPHAGAZ FLO			
<b>Explosivity or toxicity measurements</b>						
Purge or zero gas	Air, N <sub>2</sub>		ALPHAGAZ 1			
Purge or zero gas	Air, N <sub>2</sub>		ALPHAGAZ FLO			
<b>Total hydrocarbon (THC) measurement</b>						
Gas for flame	H <sub>2</sub>		ALPHAGAZ 1	ALPHAGAZ 2		
Gas for flame	H <sub>2</sub>		ALPHAGAZ FLO			
Gas for flame	H <sub>2</sub> /He		ALPHAGAZ MIX			
Gas for flame	Air		ALPHAGAZ 1	ALPHAGAZ 2		
Gas for flame	Air		ALPHAGAZ FLO			
<b>Measurement of O<sub>2</sub> using electrochemical cells coulometric or galvanic or zircon cells</b>						
Purge or zero gas	N <sub>2</sub>	N/A	ALPHAGAZ 1			
Purge or zero gas	N <sub>2</sub>	N/A	ALPHAGAZ FLO			
Zero gas	Air	N/A	ALPHAGAZ 1			
Zero gas	Air	N/A	ALPHAGAZ FLO			
<b>Measurement of O<sub>2</sub> using paramagnetic cells</b>						
Purge or zero gas	N <sub>2</sub>	ALPHAGAZ 1	N/A			
Purge or zero gas	N <sub>2</sub>	ALPHAGAZ FLO	N/A			
Zero gas	Air	ALPHAGAZ 1	N/A			
Zero gas	Air	ALPHAGAZ FLO	N/A			
<b>Measurement of H<sub>2</sub> using electrochemical cells</b>						
Purge or zero gas	N <sub>2</sub>	ALPHAGAZ 1	N/A			
Purge or zero gas	N <sub>2</sub>	ALPHAGAZ FLO	N/A			
<b>Measurement of carbon using electrochemical cells</b>						
Purge or zero gas	N <sub>2</sub>	ALPHAGAZ 1	N/A			
Purge or zero gas	N <sub>2</sub>	ALPHAGAZ FLO	N/A			

## Reference

- Gastiger M, Jurcik B. Choosing the proper gas and gas equipment for the laboratory, part 1: relationship between gas purity and detection limits in gas chromatography. *Int Lab* 1999; 29(4A):18.

*Dr. Gastiger is with the Pure Gases and Mixtures Analysis and Laboratory Market, and Dr. Jurcik is Group Manager, Application of Gases for Electronics and Laboratories, AIR LIQUIDE Corp., Claude-Delorme Research Center, BP 126, 78354 Jouy-en-Josas, Cedex, France; tel.: (33) 1 39 07 62 44; fax: (33) 1 39 07 64 66; e-mail: benjamin.jurcik@airliquide.com.*